

ORTHO MINOR SURGERY PACK




SWASTIMED
INDUSTRIES

Maximum Protection and Infection Control

www.swastimed.com


Ortho Minor Surgery Pack




Lamino Spinal
Drape




C-Arm Cover




Implant Trolley
ABS Drape




Hand Towel



Cautery Bag



Mayo ABS Cover



Instrument Trolley
ABS Drape




Surgical
Wraparound Gown




Small
Trolley Drape



Large Adhesive
Drape



Suction Bag



Side Drape



Surgical Gowns



Surgical Packs



Surgical Drapes



Accessories



**Contract
Manufacturing**

Package particulars

- 4 x Surgical Wraparound Gowns
- 1 x Lamino Spinal Drape
- 1 x Implant Trolley ABS Sheet
- 1 x Instrument Trolley ABS Sheet
- 1 x Mayo Trolley ABS Cover
- 1 x Small Trolley Sheet
- 1 x C-Arm Cover
- 2 x Large Drapes
- 2 x Side Drapes
- 1 x Suction Device Bag
- 1 x Cautery Bag
- 8 x Hand Towels

Package Information

SwastiMed's Medical Disposables are designed and manufactured with the highest level of quality and safety standards to protect the user from transfer of micro-organisms and bodily fluids while performing their duties.

All our surgical packs are designed and manufactured to meet global standards and ensure utmost safety and comfort draping at a healthcare setting. Our products are compliant with EN 13795 standard specifics.

About Us

SwastiMed Industries is a comprehensive (OEM) manufacturer of quality disposable medical supplies serving a notable clientele across the globe, protecting millions of lives.

Our state-of-the-art manufacturing unit in Chennai is well-equipped to meet bespoke requirements while conforming to highest hygiene standards designed to comfort fit, be application specific and meet safety needs.




Acquired by 30k+ healthcare professionals everyday



Trusted application in 1 million+ surgeries



Protecting more than 50 million lives



Leveraging 50+ years of combined experience



SWASTIMED
INDUSTRIES

Maximum Protection and Infection Control

99/9, Agaram Main Road, Paduvanchery, Madambakkam, Agaramthen,
Chennai - 600 126

Timings: Mon - Sat | 8 AM - 6 PM

+91 96770 88865

info@swastimed.com

www.swastimed.com